

# Archimedes' principle



Physics

Mechanics

Mechanics of liquids &amp; gases



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time

10 minutes



## Teacher information

### Application



Determination of the amount of water displaced by a mass

According to the Archimedean principle, the following statement applies:

"The static lift  $F_A$  of a body in a medium is as great as the weight of the medium displaced by the body."

Conversely, this means that a body of a certain mass  $m_K$  and a lower density than water into a container filled with water, displaces exactly the amount of water that corresponds to its weight. I.e. the weight force  $F_G$  of the body corresponds to the product of volume  $V_W$  and density  $\rho_W$  of the displaced water with the acceleration of gravity  $g$ :

$$F_G = m_K \cdot g = V_W \cdot \rho_W \cdot g \text{ [N]} \Rightarrow m_K = V_W \cdot \rho_W \text{ [kg]}$$

## Other teacher information (1/2)

**PHYWE**  
excellence in science

### Prior knowledge



Students should already have gained a basic understanding of the weight force of bodies. The students should also know the density of water.

### Scientific principle



A body immersed in a resting fluid displaces a certain amount of the fluid. This amount of fluid corresponds exactly to the weight of the immersed body if the density of the body is less than the density of the fluid. If the density is higher, the displaced volume is equal to the volume of the body and the weight force is reduced by the opposite buoyancy force, which is equal to the weight force of the displaced fluid.

## Other teacher information (2/2)

**PHYWE**  
excellence in science

### Learning objective



The students should be able to use the weight force  $F_G$  of a body first in air and then in water. From the difference they should determine the buoyancy force  $F_A$  calculate. In a second test they shall determine the mass of water displaced during immersion, calculate its weight force from this and compare it with the buoyancy force.

#### Comment:

The term "(slot)weight" is incorrect in that it refers to a mass which only receives a "weight", i.e. a weight force, under the influence of the acceleration due to gravity. The expression "piece of mass" is better.

### Tasks



#### Hint:

Stead of the beam balance, another balance with a weighing capacity of 100 g and a division of 1 g (or better) can be used.

## Safety instructions

**PHYWE**  
excellence in science

The general instructions for safe experimentation in science lessons apply to this experiment.

**PHYWE**  
excellence in science

## Student Information

## Motivation

**PHYWE**  
excellence in science



Ship floating in the sea

Due to the Archimedean principle it is possible that hot air balloons fly or ships float on the water. For this purpose, the vehicles are constructed so that the average density is lower than the medium in question. If the density of the body exceeds that of the medium, the body sinks to the ground, but its weight force is reduced by the opposite buoyancy force.

In this experiment you will learn to what extent the weight force is reduced by the buoyancy force and what the relationship is between the amount of water displaced and the weight force of the body in question.

## Tasks

**PHYWE**  
excellence in science




Does a force act on a body when immersed in water?

- Determine the weight of a body first in air and then in water.
- Determine the amount of water displaced when the body is immersed and its weight.


## Equipment

Position	Material	Item No.	Quantity
1	Weight holder, 10 g	02204-00	1
2	Slotted weight, black, 10 g	02205-01	4
3	Slotted weight, black, 50 g	02206-01	2
4	Spring balance, transparent, 2 N	03065-03	1
5	Overflow vessel 250 ml	02212-00	1
6	Beaker, 100 ml, plastic (PP)	36011-01	1
7	Graduated cylinder, 50 ml, plastic	36628-01	1
8	Pipette with rubber bulb	64701-00	1
9	Balance pan, plastic	03951-00	2
10	Lever	03960-00	1
11	Pointer for lever	03961-00	1
12	Plate with scale	03962-00	1
13	Support base, variable	02001-00	1
14	Support rod, stainless steel, l = 250 mm, d = 10 mm	02031-00	1
15	Boss head	02043-00	1
16	Holding pin	03949-00	1
17	Set of precision weights, 1g-50g	44017-01	1

## Set-up (1/4)



Connecting the tripod foot



Mounting the support rod

Build a tripod with the tripod foot, the tripod rod and the double socket.

First connect the two halves of the tripod foot to form a tripod foot.


Then screw the stand rod into the stand foot and mount the double socket.




Attaching the double socket to the stand rod

## Set-up (2/4)

**PHYWE**  
excellence in science



Fasten scale to the beam balance



Fasten the retaining bolt in double socket

Fasten the plate with scale together with the pointer using the retaining bolt in the middle of the beam scale.

Then fix the retaining bolt in the double socket.

## Set-up (3/4)

**PHYWE**  
excellence in science



Assembling the weighing pans



Taring the balance

Put the pans together and hang them on the ends of the beam balance.  
Adjust the pointer so that it points exactly to the zero mark.

## Set-up (4/4)

**PHYWE**  
excellence in science



Filling the overflow vessel with water

Fill the overflow vessel with water until it runs straight into the beaker.


Wait until no more water drips out, then dry the beaker carefully.



## Procedure (1/4)

**PHYWE**  
excellence in science

- Determine the mass  $m_0$  of the dry beaker with the beam balance and note the value in the protocol.



Balanced beam scales:

Determination of the mass of the empty cup

## Procedure (2/4)

**PHYWE**  
excellence in science



Determination of the weight forces in air

- Determine the weight in air for the masses 50 g, 100 g and 150 g with the dynamometer  $F_{G_{iL}}$ .
- Note the measured values in the table in the protocol.

## Procedure (3/4)



Determination of the weight forces in water

- Place the well-dried beaker under the maximum-filled overflow vessel and immerse the weight plate with the pieces of mass for the total masses of 50 g, 100 g and 150 g completely in the overflow vessel one after the other.
- Read the value for the weight in water  $F_{G,W}$  for the masses.
- Wait until no more water drips out and then determine the mass of the displaced water including the beaker  $m_1$  with the beam balance.
- Enter all measured values in the table.
- Note: Make sure for each total mass that the beaker is dry and the overflow vessel is filled to the maximum.

## Procedure (4/4)

- To disassemble the tripod base, press the buttons in the middle and pull both halves apart.



Disassembling the tripod base



# Report

## Table

Enter them in the appropriate fields.

Determine the buoyancy forces of the masses:  $F_A = F_{G_{iL}} - F_{G_{iW}}$

$m_0 =$

Determine the mass of the overflowed water:  $m_W = m_1 - m_0$

Calculate the associated weight forces from this:  $F_W = m_W \cdot g (g = 9,81m/s^2)$

m [g]	$F_{G_{iL}}$ [N]	$F_{G_{iW}}$ [N]	$F_A$ [N]	$m_1$ [g]	$m_W$ [g]	$F_W$ [N]
50						
100						
150						

## Task 1

Compare the results for  $F_A$  with those of  $F_W$ . What do you find?

- $F_W$  is greater than  $F_A$ .
- $F_A$  is greater than  $F_W$ .
- $F_A$  and  $F_W$  are the same size.

✓ Check

## Task 2

With which two methods can the buoyancy force  $F_A$  determine?

- The buoyancy force can be determined directly by measuring the mass in water in terms of its weight.
- One can determine the buoyancy force by building the difference of the weight forces in air and water.
- One can determine the buoyancy force by determining the weight of the displaced water.

✓ Check

### Task 3

**PHYWE**  
excellence in science

How does the buoyancy force affect a submerged body?

- It counteracts its weight and thus makes it appear lighter.
- It has no effect on the body.
- It works together with its weight and makes it appear heavier.

✓ Check

### Task 4

**PHYWE**  
excellence in science


When does a body swim, when does it sink?

- A body floats when it is gently laid on the water.
- A body always floats when its average density is lower than that of water.
- A body always floats when the effective buoyancy force is greater than its weight in water

✓ Check

Slide	Score/Total
Slide 20: settlement of $\backslash(F_A)$ and $\backslash(F_W)$	0/1
Slide 21: Determination of the buoyancy force	0/2
Slide 22: impact of $\backslash(F_A)$	0/1
Slide 23: buoyancy	0/2

Total amount  0/6

 Solutions

 Repeat

 Exporting text