

Image obtained with a convex lens

Task and equipment


Information for teachers

Image obtained with a convex lens

Task and equipment

Task

Equipment



Position No.	Material	Order No.	Quantity
1	Support base, variable	02001-00	1
2	Object -L-, glass bead	11609-00	1
3	Lens on slide mount, $f=+100\text{mm}$	09820-02	1
4	Slide mount for optical bench	09822-00	1
5	Screen, white, 150x150mm	09826-00	1
6	Support rod, stainless steel, $l = 600 \text{ mm}$, $d = 10 \text{ mm}$	02037-00	2
6	Meter scale for optical bench	09800-00	1
7	Bottom with stem for light box	09802-10	1
8	Light box, halogen 12V/20 W	09801-00	1
	PHYWE power supply DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1

Set-up and procedure

Set-up

Using the two support rods and the variable support base assemble the optical bench and place the scale for the optical bench against the front support rod.



Fig. 1



Fig. 2

Place the bottom with stem under the light box.



Fig. 3



Fig. 4

Clamp it in the left part of the support base so that the lens end points away from the optical bench.



Fig. 5

Insert an opaque cover in front of the lens and the cover with the glass-bead-L into the slot at the other side of the light box.



Fig. 6

Complete the setup by mounting the lens and screen according to the figure 7.




Fig. 7

Procedure

Connect the light box to the power supply (12 V AC) and switch it on.

Student's Sheet

Printed: 18.04.2017 09:54:54 | P1068200



Fig. 8

Set the distance between the lens and the glass-bead-L (the object distance g) so that $g > 2f$. Then move the screen until the -L- is sharply focused.




Fig. 9

Measure the image distance b and compare it with the focal length f ; observe the image. Record your values in table 1 in the report. To characterise the three major properties of the image use the words: "upright" or "inverted"; "enlarged", "reduced" or "same size"; "real" or "virtual".




Fig. 10

Perform the same steps for the other cases listed on the task page and note your results in table 1 in the report. Switch off the power supply.

Report: Image obtained with a convex lens

Result - Table 1

Note down your observations in table.

Object distance	Image distance	Properties of the image		
$g > 2f$	$2f > b > f$	1	1	1
$g = 2f$	$\frac{1}{\pm 0}$	inverted	1	1
$2f > g > f$	$\frac{1}{\pm 0}$	1	enlarged	1
$g < f$	negative	1	1	virtual

Evaluation - Question 1

Under what conditions is the real image formed by a convex lens?

.....

.....

.....

.....

Evaluation - Question 2

Under what conditions is a virtual image always formed by a convex lens?

.....

.....

.....

.....

Evaluation - Question 3

What happens when the object is located at the focal length of the convex lens ($g = f$)?

.....

.....

.....

.....

Evaluation - Question 4

What can you say about the images formed by a convex lens which a piece (for example, the lower haft) is the broken off?

.....

.....

.....

.....

Evaluation - Question 5

There is a simple optical device which one often uses when one desires to examine small objects or see details of a large object. Surely you have used this device.

- What is it called?
- Give several examples for its use.

.....

.....

.....

.....