


# Properties of solar cells - the dependence on the illuminating intensity (Item No.: P1382600)

## **Curricular Relevance**




## **Principle and equipment**

## Principle

The dependence of the idle voltage generated by a solar cell and also of the short-circuit current on the illuminating intensity is to be examined.





**PHYWE** excellence in science

Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 info@phywe.de www.phywe.com

## Equipment

| Position No. | Material                                | Order No. | Quantity |
|--------------|---|-----------|----------|
| 1            | Multimeter ADM2, demo., analogue        | 13820-01  | 1        |
| 2            | Demo Physics board with stand           | 02150-00  | 1        |
| 3            | Ceramic lamp socket E27                 | 06751-01  | 1        |
| 4            | Solar cell (2.5x5)cm,module DB          | 09470-00  | 1        |
| 5            | Clamp on holder                         | 02164-00  | 1        |
| 6            | Connector interrupted, module DB        | 09401-04  | 1        |
| 7            | Electr.symbols f.demo-board,12pcs       | 02154-03  | 1        |
| 8            | Connector, straight, module DB          | 09401-01  | 2        |
| 9            | Connector, angled, module DB            | 09401-02  | 4        |
| 10           | Connector, T-shaped, module DB          | 09401-03  | 2        |
| 11           | Filament lamp,220V/120W,w.refl.         | 06759-93  | 1        |
| 12           | Scale for demonstration board           | 02153-00  | 1        |
| 13           | Bosshead, turnable                      | 02048-04  | 1        |
| 14           | Support rod, stainless steel, 500 mm    | 02032-00  | 1        |
| 15           | Connecting cord, 32 A, 1000 mm, red     | 07363-01  | 1        |
| 16           | Connecting cord, 32 A, 1000 mm,<br>blue | 07363-04  | 1        |
|              | Preliminary Experiment:                 |           |          |
|              | Motor, 2 V-                             | 11031-00  | 1        |
|              | Disc for motor, 2 V-                    | 11031-01  | 1        |
|              | Clamp on fixing magnet                  | 021541-01 | 1        |

**PHYWE** excellence in science

Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 info@phywe.de www.phywe.com Printed: 09/11/2017 08:43:08 | P1382600

## advanced

## Set-up and procedure

The solar cell is illuminated by a reflector lamp. The lamp is moved to the side and is realigned to face the solar cell at each position. The intensity of the illumination is thereby changed. The displacing of the lamp represents the daily movement of the sun. The idle voltage and short-circuit current are measured.

#### **Qualitative preliminary experiment**


The performance of a solar cell is dependent on the intensity of the illumination. This can be demonstrated with the help of a small motor.

- Fix the motor to the magnetic board with the clamp, attach the disc.
- Connect the motor to the solar cell.
- Change the illuminating intensity, starting with the greatest distance away of the lamp, i.e. with the smallest illuminating intensity.
- Observe the motor

The motor only runs slowly under weak illumination. The stronger the illumination, the quicker it runs. This is particularly true for small illuminating intensities, so that these should be used to start with. A further point to consider is that the motor keeps running for some time, even when the lamp has already been switched off, for example.

#### Measurement of the idle voltage

- Connect up the circuit as in Fig. 1; place the modules at the left edge of the board and fix the lamp at the top edge of the board, vertically above the solar cell, using the support material and the clamp; bring the zero mark on the scale to the position of the lamp.
- Select the 1 V- measurement range, measure the idle voltage  $U_0$ , enter the value in Table 1.
- Move the lamp to the side in steps of 10 cm and at each distance b realign the lamp to face the solar cell; note the idle voltage for each position.



#### Measurement of the short-circuit current

- Change the circuit to that shown in Fig. 2; select the 300 mA- measurement range, measure the short-circuit current  $I_{\rm SC}$ , and note the measured value.
- Reduce the distance b of the lamp from the solar cell in steps of 10 cm; realign the lamp at each distance, measure and note the short-circuit current.




Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

### **Student's Sheet**

Printed: 09/11/2017 08:43:08 | P1382600





## **Observation and evaluation**

## Observation

| Table 1        |                  |                  |
|----------------|------------------|------------------|
| $\frac{b}{cm}$ | $\frac{U_0}{cm}$ | $\frac{I_K}{mA}$ |
| 0              | 0.56             | 270              |
| 10             | 0.54             | 250              |
| 20             | 0.54             | 200              |
| 30             | 0.52             | 170              |
| 40             | 0.51             | 115              |
| 50             | 0.50             | 92               |
| 60             | 0.49             | 70               |
| 70             | 0.48             | 48               |
| 77             | 0.48             | 36               |

## **Evaluation**

The graphs shown in Fig. 3 and Fig. 4 are given on plotting the measured values. The following can be derived from them:

The idle voltage of a solar cell only decreases slightly with a decrease in the illuminating intensity. In contrast to this, a decrease in the illuminating intensity has a considerable effect on the short-circuit current of a solar cell. The idle voltage  $U_0$  is the maximum voltage generated by a current source when it is under no load; the short-circuit current is the maximum value of the current of a current source.

