


Chromoplasts



Biology

Microscopy / Cell Biology

Basics of Microscopy & Work Technology

Biology

Microscopy / Cell Biology

Plants & Fungi

Biology

Microscopy / Cell Biology

Cell structure

Nature & technology

From the very small & the very big

Nature & technology

Plants & animals



Difficulty level

easy



Group size

1



Preparation time

10 minutes




Execution time

30 minutes



Teacher information

Application



Cells of peppers at low magnification

The fruits and blossoms of many plants have strong, bright colours. The coloured fruits attract animals, are eaten by them and the seeds are excreted elsewhere. This is how the plant is spread. Coloured flowers attract insects, which harvest the nectar from the flowers. In addition, the pollen is transported from plant to plant, thus ensuring fertilization. Yellow and red dyes are usually found in certain cell organelles, the chromoplasts.

Other teacher information (1/3)

PHYWE
excellence in science

Prior knowledge



The cytoplasm of plant cells contains numerous plastids surrounded by thin membranes. They are usually ovoid in shape and are distinguished by their colour. Due to photosynthesis, the most important group is that of chloroplasts. The chromoplasts contain carotenes or xanthophylls. Some of the colours have a meaning for attracting animals (e.g. flowers - insects, fruits - birds).

Scientific principle



The students should understand the preparation instructions as a stimulus. They can make peel preparations, cut preparations, but also scraping preparations of carrots or squeezing preparations of the soft flesh of tomatoes. It is important to make sure that the preparations become thin, i.e. that they are translucent.

Other teacher information (2/3)

PHYWE
excellence in science

Learning objective



The students learn how chromoplasts are present in form and arrangement in the cell.



Tasks



Making and microscoping the preparation


1. Red pepper - fruit skin
2. Capuchin Cress - Petal
3. Rose - pulp

Other teacher information (3/3)



Chromoplasts are present in all yellow, red and orange parts of the plant. The preparations suggested on the student worksheet are particularly suitable, but students should be encouraged to examine fruits and flowers that are easily accessible to them. The pupils can bring coloured vegetables from the kitchen at home and flowers from the garden. Blue and purple flowers are not suitable, as the most common dye anthocyanin is dissolved in the vacuoles. Fruits of the rose can be harvested from wild roses or from bed roses.

Safety instructions




- Razor blades are sharp. Extreme caution is advised when handling them.
- Working with microscopes for too long can lead to physical discomfort (fatigue, headaches, nausea), especially when students are inexperienced.
- Microscopes are sensitive. During transport and handling, care should be taken to ensure that everything is done carefully and without rushing.
- The general instructions for safe experimentation in science teaching apply to this experiment.



Student Information


Motivation




Cells of peppers at low magnification

In this experiment you will learn how to recognize the chromoplasts in different preparations of different plants and fruits and how to describe their shape and arrangement.

Tasks



Fruits of the rose (400x)



Peppers (400x)

Preparation and microscopy of preparations

1. Red pepper - fruit skin
2. Capuchin Cress - Petal
3. Rose - pulp

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Binocular student microscope, 1000x, mechanical stage	MIC-129A	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs	64685-00	1
4	Beaker, 100 ml, plastic (PP)	36011-01	1
5	Dropping pipette with bulb, 10pcs	47131-01	1
6	Tweezers, straight, pointed, 120mm	64607-00	1
7	Scalpel holder	64615-00	1
8	Scalpel blades, rounded tip, 10 off	64615-02	1

Procedure (1/3)

PHYWE
excellence in science




(1) Prepare the sample: Red pepper - fruit skin

- Prepare a slide with a drop of water.
- With the tweezers a piece of the skin of the pepper fruit is peeled off.
- The cuticle is placed directly into the water drop, covered and microscopically examined.

Procedure (2/3)

(2) Prepare the sample: Capuchin cress - petal

- Prepare a slide with a drop of water.
- A yellow or red petal is placed over the finger and a thin surface incision is made with the scalpel.
- The preparation is placed directly into the water drop, covered and microscopically examined.



Procedure (3/3)

PHYWE
excellence in science

(3) Prepare the sample: Rose - pulp

- Prepare a slide with a drop of water.
- From the red fruit of a rose (rosehip) a very thin piece of fruit is cut off. The even flesh of the fruit is examined, not the seed.
- The preparation is placed directly into the water drop, covered with a cover glass and microscopically examined.

PHYWE
excellence in science

Report

Task 1

PHYWE
excellence in science

In the plant parts with which colours are chromoplasts contained?

 Greens Oranges Yellow Red Check

Task 2

PHYWE
excellence in science


What is a possible meaning of the colours in the chromoplasts?

 Defence Attracting animals (e.g. flowers - insects; fruits - birds) Sex attractant Lights at night

Task 3

PHYWE
excellence in science


Microscopy with increasing magnification. Drawing from the specimen **Red pepper - fruit skin** a cell. Use a coloured pencil to show the location and number of chromoplasts.



Task 4


PHYWE
excellence in science

Microscopy with increasing magnification. Drawing from the specimen **Capuchin Cress - Petal** a cell. Use a coloured pencil to show the location and number of chromoplasts.



Task 5

Microscopy with increasing magnification. Drawing from the specimen **Rose - pulp** a cell. Use a coloured pencil to show the location and number of chromoplasts.



Slide

Score/Total


Slide 15: Plant parts with chromoplasts

0/3

Slide 16: Colours in chromoplastics

0/1

Total amount

 Solutions Repeat