Fish Gills



Biology	Microscopy / Cell Biology			
Biology	Microscopy / Cell Biology		Humans & Animals	
Biology	Animal Physiolog	gy / Zoology		
Nature & technology		From the very small & t	he very big	
Nature & technology		Plants & animals		
Difficulty level easy	QQ Group size	Preparation 10 minu	n time tes	Execution time 30 minutes







Teacher information

Application





Many simply built animals, such as the slipper animal or the earthworm, absorb oxygen for all life processes over the entire surface. In fish, many skin folds, the so-called gills, have formed in the course of evolution. The water flows over this contact surface and the oxygen contained in it migrates through the gill skin into the blood of the fish. The fish breathes.

Two gill leaflets 40x





Other teacher information (2/4)



Learning

The students should recognize the structure of the gills, the gill arch and the gill leaflet and thus understand the function of these body parts. In doing so, they will recognize that the structure of the fish gills results in an increase in surface area.

Tasks

Students will dissect and microscope gills, gill arches, and gill leaflets.



Other teacher information (3/4)

Notes on material procurement

Fresh fish can be obtained from an angler or at the fish counter in the supermarket. In principle, all fish are suitable for preparation if they have a minimum size of approx. 10 cm so that they can also be examined macroscopically. If you have anglers in the class, they can prepare a frozen stock.



Gill leaflets 100x



Gill leaflets 400x

Other teacher information (4/4)

Notes on implementation

expose gills: Here the students have to cut vigorously for the first time, which may cost overcoming. The gill clasps consist of a bony substance.

The gill archThe pupils can recognise all the components of the bow on the white background. Depending on the didactic orientation, the parts gill trap and gill arch can also be named and explained. By drawing the macroscopic structure, an idea of the surface enlargement is introduced.

Microscopy of the gill leaflets: Students should articulate the different observations in the ascending magnifications (including magnification with a magnifying glass). The main cognitive learning objective of this lesson: students will recognize the principle of surface area magnification. More contact or reaction surface area can produce higher performance. This principle can be transferred to other areas (alveoli, chloroplasts).



Safety inst	tructions	PHYWE excellence in science
	 Working with microscopes for too long can lead to physical discomfort (feedback) 	atigue,
	headache, nausea), especially when students are untrained.	avoid
	accidents!	
	 Microscopes are sensitive. During transport and handling, care should be ensure that everything is done carefully and without rushing. 	e taken to
	 The general instructions for safe experimentation in science lessons applexperiment. 	ly to this

Student Information



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

Motivation





Two gill leaflets 40x

Many simply built animals, such as the slipper animal or the earthworm, absorb oxygen for all life processes over the entire surface. In fish, many skin folds, the so-called gills, have formed in the course of evolution. The water flows over this contact surface and the oxygen contained in it migrates through the gill skin into the blood of the fish. The fish breathes.

Tasks





Prepare gills, gill arches and gill leaflets. Explore and, if necessary, microscope the preparations. Do you recognise the principle of surface enlargement?



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Binocular student microscope, 1000x, mechanical stage	MIC-129A	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs	64685-00	1
4	Scissors,straight,pointed,I 110mm	64623-00	1
5	Magnifier, plastic, 5x, d=35mm	88002-01	1
6	Tweezers,straight,pointed,120mm	64607-00	1

Procedure (1/2)

PHYWE excellence in science

Expose the gills

• The gills are protected from the outside by the gill cover. Remove the lid with the dissecting scissors. Now you can see the gill arches.

The gill arch

• Prepare a complete gill sheet. Place the sheet on a white surface and examine it closely with the magnifying glass.



Procedure (2/2)

Microscope a gill leaflet...

- At the gill arches there are many thin membranes, the gill leaflets. Try to recognize the blood vessels under the microscope!
- Cut off a small piece of the gill skin with the dissecting scissors and place it in the water drop.
- Microscope gradually up to the highest magnification and draw in detail the gill leaflets in the record.



PHYWE excellence in science









Report

Task 1	PHYW excellence in scier
Choose the correct answers.	
Fish absorb water through their gills and release it through their mouths.	
Fish take in water through their mouths and release it through their gills.	
The structure of the gills causes an enlargement of the surface, whereby more oxygen can be taken up.	
The structure of the gills causes a surface reduction, whereby more oxygen can be taken up.	
Check	



Task 2	PHYWE excellence in science
Independent of the service of the servic	In the course of evolution, fish have formed many skin folds, the so-called gills. The water flows over this contact surface and the oxygen contained in it travels through the gill skin into the fish's blood. The fish breathes.

Task 3

Draw in detail gill arches and gill leaflets and label them.





PHYWE excellence in science

Slide	Score/Total
Slide 15: Fish	0/2
Slide 16: Multiple tasks	0/2
	Total 0/4
	Solutions Repeat