

# Franck-Hertz experiment with a Ne-tube



P2510315

Physics

Modern Physics

Quantum physics

Chemistry

Physical chemistry

Atomic structures &amp; properties



Difficulty level

easy



Group size

1



Preparation time

10 minutes




Execution time

10 minutes

**PHYWE**  
excellence in science

## General information

## Application

**PHYWE**  
excellence in science

Franck-Hertz experiment with neon gas

The Franck-Hertz experiment shows the absorption of kinetic energy of electrons by neon atoms, or by mercury atoms in the original experiment.

## Other information (1/2)

**PHYWE**  
excellence in science

### Prior knowledge



The prior knowledge required for this experiment is found in the theory section.

### Scientific principle



Electrons are accelerated in a tube filled with neon gas. The excitation energy of neon is determined from the distance between the equidistant minima of the electron current in a variable opposing electric field.

## Other information (2/2)

**PHYWE**  
excellence in science

### Learning objective



Understanding energy quanta through Franck-Hertz experiment with neon gas

### Tasks



1. Record the countercurrent strength  $I$  in a Franck-Hertz tube as a function of the anode voltage  $U$ .
2. Determine the excitation energy  $E$  from the positions of the current strength minima or maxima by difference formation.

## Safety instructions

**PHYWE**  
excellence in science

For this experiment the general instructions for safe experimentation in science lessons apply.

For H- and P-phrases please consult the safety data sheet of the respective chemical.

## Theory (1/5)

**PHYWE**  
excellence in science

Niels Bohr introduced the planetary model of the atom in 1913: An isolated atom consists of a positively charged nucleus about which electrons are distributed in successive orbits. He also postulated that only those orbits occur for which the angular momentum of the electron is an integral multiple of  $h/2\pi$ , i.e.  $n \cdot h/2\pi$ , where  $n$  is an integer and  $h$  is Planck's constant.

Bohr's picture of electrons in discrete states with transitions among those states producing radiation whose frequency is determined by the energy differences between states can be derived from the quantum mechanics which re-placed classical mechanics when dealing with structures as small as atoms.


It seems reasonable from the Bohr model that just as electrons may make transitions down from allowed higher energy states to lower ones, they may be excited up into higher energy states by absorbing precisely the amount of energy representing difference between the lower and higher states.

## Theory (2/5)

**PHYWE**  
excellence in science

James Franck and Gustav Hertz showed that this was, indeed, the case in a series of experiments reported in 1913, the same year that Bohr presented his model. Franck and Hertz used a beam of accelerated electrons to measure the energy required to lift electrons in the ground state of a gas of mercury atoms to the first excited state.

In present experiment, a tube filled with neon gas is used. The electrons emitted by a thermionic cathode are accelerated between cathode C and anode A in the tube filled with neon gas and are scattered by elastic collision with neon atoms.



Principle of the measurement.

## Theory (3/5)

**PHYWE**  
excellence in science

From an anode voltage  $U_1$  of 16,8 V, however, the kinetic energy of the electrons is sufficient to bring the valence electron of the neon to the first excitation level by an inelastic collision. Because of the accompanying loss of energy, the electron can now no longer traverse the opposing field between anode A and counter electrode S: the current  $I$  is at a minimum.

If we now increase the anode voltage further, the kinetic energy of the electron is again sufficient to surmount the opposing field: the current strength  $I$  increases.

When  $U_1 = 2 \cdot 16.8 \text{ V}$  the kinetic energy is so high that two atoms in succession can be excited by the same electron: we obtain a second minimum.

These minima are not, however, very well-defined because of the initial thermal distribution of the electron velocities.

## Theory (4/5)

The voltage  $U_1$  between anode and cathode is represented by

$$U_1 = U + (\Phi_A - \Phi_C)$$

where  $U$  is the applied voltage, and  $\Phi$  the work function voltages of the anode and cathode re-spectively. As the excitation energy  $E$  is determined from the voltage differences at the minima, the work function voltages are of no significance here.

According to the classical theory the energy levels to which the mercury atoms are excited could be random. According to the quantum theory, however, a definite energy level must suddenly be assigned to the atom in an elementary process.

The course of the  $I/U_A$  curve was first explained on the basis of this view and thus represents a confirmation of the quantum theory. The excited neon atom again releases the energy it has absorbed, with the emission of a photon.

## Theory (5/5)

When the excitation energy  $E$  is 16.8 eV, the wavelength of this photon is

$$\lambda = \frac{ch}{E}$$

where  $c = 2.9979 \cdot 10^8 \frac{m}{s}$  and  $h = 4.136 \cdot 10^{-15} eVs$

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Franck-Hertz Ne-tube w. housing</a>	09105-40	1
2	<a href="#">Franck-Hertz control unit</a>	09105-99	1
3	<a href="#">Connect.cord f.Franck-H. Ne-tube</a>	09105-50	1
4	<a href="#">Software Measure Franck-Hertz experiment</a>	14522-61	1
5	<a href="#">Converter USB - RS232, active</a>	14602-10	1
6	<a href="#">Data cable, plug/ socket, 9 pole</a>	14602-00	1
7	<a href="#">Screened cable, BNC, l = 750 mm</a>	07542-11	1

## Additional equipment

**PHYWE**  
excellence in science

Position	Material	Quantity
1	PC	1

**PHYWE**  
excellence in science

## Setup and procedure



## Setup


 PHYWE  
 excellence in science


Experimental setup

Set up the experiment as shown in the figure. For details see the operating instructions of the unit 09105-99.


Connect the Franck-Hertz operating unit to the computer port COM1, COM2 or to USB port (use USB to RS232 Adapter Converter 14602-10).

## Procedure


 PHYWE  
 excellence in science


Start the measure program and select Franck- Hertz experiment Gauge. The window "Frank- Hertz-experiment – measuring" appears. The optimum parameters are different for each Ne-tube. You find the specific parameters for your device on a sheet which is enclosed in the package of the Ne-tube.

Choose the parameters for U1, U2, U3 and UH as given on that sheet and make sure that the rest is set as shown in Fig. 2. Press the continue button.



Measuring parameters.

## Evaluation (1/2)



The graph of  $I/U_1$  shows equidistant maxima and minima.

For our evaluation we determine the voltage values of the minima. From the differences between these values we obtain the excitation energy  $E$  of the neon atom by taking an average. By evaluating the measurements in the figure we obtained the value

$$E = (17.4 \pm 0.7) \text{ eV}$$

Example of a Franck-Hertz curve recorded with Ne-tube.

## Evaluation (2/2)

Fill in the blanks:

In a tube filled with \_\_\_\_\_, electrons are emitted by glow emission from a thermionic \_\_\_\_\_ and accelerated to the \_\_\_\_\_ by the acceleration voltage. If the voltage is increased slowly, the measured current values initially \_\_\_\_\_ exponentially, up to a certain voltage and then \_\_\_\_\_ slowly. The current values \_\_\_\_\_ again when voltage is twice the value of the initial voltage (at which the current rise for the first time).

- anode
- increase
- cathode
- decrease
- increase
- neon gas

✓ Check

Slide

Score / Total

Slide 17: Summary experiment

0/6

Total Score



 Show solutions

 Retry