## **Dependence of the reaction behavior of metals**



Chemistry	Inorganic chemistry	Chemistry	of metals
Difficulty level	<b>QQ</b> Group size	C Preparation time	Execution time
easy	2	10 minutes	10 minutes








# **Teacher information**

## **Application**







**Test setup** 


As with other substances, metals react at different rates when heated in air. The speed and intensity of the reaction depends on the one hand on the type of metal (e.g. noble or base metal), and on the other hand on its degree of fragmentation. With a higher degree of fragmentation, the surface area of the substance becomes larger and the probability of an "encounter" with the reactant increases. Substances with a high degree of fragmentation have a larger reactive surface. This means that they react faster and more violently than substances with a lower degree of fragmentation. The violence and speed of the reaction also depends on the type of metal. Here the terms "noble" and "base" play a major role. Base metals react well. In contrast, noble metals react less well or not at all when heated in air.



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107











# **Student Information**

Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

## Motivation






Candy sugar has a lower degree of fragmentation

The degree of fragmentation influences the speed and violence of a reaction. The degree of fragmentation is the ratio of the surface area to the volume of a substance. To better understand the term, the following examples from everyday life.

Wood chips are easier to light than a thick log, iron powder easier than iron nail.

Another example of the dependence of the reaction rate on the degree of disintegration is the decomposition of foodstuffs. Candy sugar has a lower degree of disintegration compared to fine crystal sugar. Candy sugar (pieces) therefore take longer to dissolve in water, even with the same quantities.






#### Equipment

Positio	on Equipment	ltem no.	Quantity															
1	Evaporating dish, 75 ml, 3 top-d = 80 mm	32516- 00	1	Wire mesh 2 with ceramic, 160 x 160 mm	01	Burning 13 spoon (phosphorus spoon)	33346- 00	Crucible tongs, 1 4 stainless steel, I = 200 mm	33600-	Watch glass 15 bowl, d = 60 mm	84570- 2 00 2	Safety glasses "classic" - OneSize, Unisex	Powder spatula, 17 steel, I = 150 mm	47560- 00	Iron, 18 coarse powder, 500 g	30067- 50	Platinum wire, d = 190.3 mm, I = 100 mm	

### Set-up






- Take a gas burner and connect the gas burner to the gas supply.
- Follow your teacher's instructions when connecting the gas burner.
- $\,\circ\,\,$  Place the burner in the center of the workstation.
- Now take two watch glass bowls.
- Put a little metal powder (a spatula tip) on each watch glass.
- Mark the watch glass shells with a laboratory recorder.

## Procedure (1/2)



Put a spatula tip of iron powder into the combustion spoon (Fig. bottom left).


Set the non-luminous burner flame and heat the iron powder vigorously. Hold the spoon at a slight angle so that no powder falls out.





PHYWE excellence in science

## Procedure (2/2)



- In between, let the burn spoon cool down and clean it from metal traces.
- $\circ\;$  Take some iron powder from the watch glass dish with the spatula.
- Hold the burner at an angle with the left hand and tap the spatula lightly above the flame from a sufficient height so that some iron powder enters the flame.
- Use only the smallest amounts of powder and make sure that it does not get into the burner opening.
- Grasp the platinum wire with the crucible tongs and heat it intensively in the hottest part of the burner flame.
- $\circ~$  Look closely at the wire before and after the test.

**Disposal**: Add the metal powder in the evaporating dish to the heavy metal waste.





## Report




Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

Task 1			
	Write down	your observations!	
	Metal	Burning spoon open flame	
	Iron powder		
7	Zinc powder		
	Platinum wi		
		· · · · · · · · · · · · · · · · · · ·	

Task 2			<b>PHYWE</b> excellence in science
Complete the sentences.			
The metal powders react much	violently in the open f	flame than when heated in	the same
the combustion spoon. They partially burn in the	ne open flame under	heat	burn
development, which can be seen from the brig	ht flame.		different
Sheets, on the other hand, do not	and react very	compared to the	more
powder. However, the substances formed in th	e reactions are	, even though the	slowly
reaction rates were			strong
Check			

Task 3	<b>PHYWE</b> excellence in science
Complete the sentences.	
1. the platinum wire change even after long and strong heating. So no chemical	noble
reaction takes place here. Platinum is a metal that does react when heated in air.	fine does not
2. The reactivity of metals when heated in air depends on the type of metal. metals	not noble
react well, metals react less well or not at all. Like metals react better the they are distributed.	precious
Check	





Slide				Score/Total				
Slide 15: Metal pow	der combustion			0/6				
Slide 16: Reactivity of metals				0/6				
Slide 17: Burn reac	tion rate			0/1				
			Total	0/13				
	Solutions	2 Repeat	Export text					