
Hardness of water





PHYWE excellence in science



Teacher information


Application




Water for body cleansing

Water is not only a necessary foodstuff (drinking water), but also plays in everyday life a special position. (Drinking) water forms the basis of almost all beverages. We also use water to wash and clean ourselves and objects. Water is also used as a solvent in many technical processes. In some of these processes (e.g. in the washing machine), the "quality" of the water plays an essential role. This can be seen in everyday life by the calcification of dishes or kettles. The degree of calcification is directly linked to the so-called water hardness. In principle, the higher the water hardness of the water, the more calcified it is, the quicker it calcifies machines or everyday products such as tableware.





Other teacher information (2/2)



excellence in science

Safety instructions

PHYWE excellence in science

The general instructions for safe experimentation in science lessons apply to this experiment.

For H- and P-types please consult the safety data sheet of the respective chemical.





Student Information



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

Motivation



PHYWE excellence in science



Water is not equal to water, so we divide water into hard and soft water. The water hardness is determined by magnesium and calcium ions in Water. Basically: The higher the concentration of these metal salts, the harder the water. The negative consequences of a high water hardnessärte we know from everyday life. High water hardness is responsible for the fact that appliances and dishes calcify. In addition, a high water hardness reduces the washing effect of detergents. In order to give a more scientific definition of water hardness, the water hardness was defined in a temporary and permanent Hardness classified. The temporary hardness refers to sich only refers to the proportion of carbonate ions in Water.





1. In this experiment you examine different water samples for their (different) water hardness and how the water changes according to the proportion of (above all) calcium and magnesium salts into hard und soft water to divide.

- 2. Experimentally prove that lime (scale) precipitates when hard water is heated.
- 3. Check whether the so-called lime soap (which has precipitated out of the solution) can be dissolved by adding water softener.



Equipment

| Position | Material | Item No. | Quantity |
|----------|---|----------|----------|
| 1 | Support base, variable | 02001-00 | 1 |
| 2 | Support rod, stainless steel, I=370 mm, d=10 mm | 02059-00 | 1 |
| 3 | Wire gauze with ceramic, 160 x 160 mm | 33287-01 | 1 |
| 4 | Beaker, Borosilicate, low form, 250 ml | 46054-00 | 1 |
| 5 | Beaker, Borosilicate, tall form, 250 ml | 46027-00 | 1 |
| 6 | Ring with boss head, i. d. = 10 cm | 37701-01 | 1 |
| 7 | Protecting glasses, clear glass | 39316-00 | 1 |
| 8 | Glass rod, boro 3.3, I=200mm, d=5mm | 40485-03 | 1 |
| 9 | Spatula, powder, steel, I=150mm | 47560-00 | 1 |
| 10 | Butane burner with cartridge, 220 g | 32180-00 | 1 |
| 11 | Soap solu.(Boutron-Boudet) 250 ml | 30221-25 | 1 |



PHYWE excellence in science

Additional equipment

Position Equipment Quantity

1

- 1 Water samples 1
- 2 Water softener



- Assemble the tripod from the tripod base and the tripod rod. Note the two upper pictures in the illustration on the right
- Attach the tripod ring to the tripod rod and place the wire netting on it. Please note the lower two illustrations in the picture on the right.











PHYWE excellence in science

PHYWE excellence in science

PHYWE excellence in science

Set-up (2/2)

- Take two beakers and label them with 1 and 2.
- Fill one beaker 1 two-thirds full with tap water, the other two-thirds full with distilled water (as shown in the figure on the right)



Beakers with different water samples

Procedure (1/2)



- Pour about 5 ml of soap solution into both beakers (see illustration top left).
- Stir with the glass stick.



- $\circ\;$ Leave the beakers to stand still for a short time after stirring.
- $\circ~$ Note down what was observed during stirring.



Procedure (2/2)

Place the beaker with the tap water on the wire netting (bottom left figure) and heat it to boiling point (bottom middle figure). Adjust the burner flame so that the water just boils. Add half a spoonful of water softener to the boiling water (bottom right illustration)







Report



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107

| Task 1 | PHYWE excellence in science |
|---|---------------------------------------|
| Describe your observations | |
| When stirring after adding soap solution, foam is produced in both beakers. However, foam formation is much stronger when using than with | water soltener |
| In the beaker containing tap water, are formed after a certain time and settle. After adding the these | insoluble flakes |
| flakes dissolve again and the water stops boiling. | distilled water |
| ✓ Check | |
| | |

Task 2

| Since tap water does not foam as much as distilled |
|--|
| water, the former must contain dissolved |
| substances that react with the soap. How can this be |
| seen? |

| 0 | The tap water starts boiling | |
|---|------------------------------|--|
| | | |

- $\ensuremath{{\ensuremath{\mathsf{O}}}}$ At the formation of the precipitate (lime soap)
- O Soap foams equally in both solutions





Water samples in comparison



PHYWE excellence in science

| ide | Score/Total |
|--------------------------------------|------------------|
| lide 16: Water hardness | 0/4 |
| lide 17: Water samples in comparison | 0/ 1 |
| | Total amount 0/5 |
| | |
| | |
| | |
| | |
| | |
| Solutions | C Repeat |