The swelling







Teacher information

Application



Experimental setup for pea swelling



Swelling is generally defined as the absorption of liquid or vapour by a swelling body with an increase in its volume. Swelling is a purely physical process, which may also be observed, for example, when gelatin is soaked in water. However, even air-dried seeds, whose water content is usually around 5%, swell when they come into contact with water. Again, this is a purely physical aspect in which the metabolism of the seeds is not directly involved. For example, seeds which are already dead and no longer germinable swell just as well as those which are germinable.



Robert-Bosch-Breite 10 37079 Göttingen





Robert-Bosch-Breite 10 37079 Göttingen

www.phywe.de

Safety instructions





• The general instructions for safe experimentation in science lessons to be applied to this experiment.





Student Information



Tel.: 0551 604 - 0 Fax: 0551 604 - 107

Motivation





Swelling is generally defined as the absorption of liquid or vapour by a swelling body with an increase in its volume. Swelling is a purely physical process, which may also be observed, for example, when gelatin is soaked in water. However, even air-dried seeds, whose water content is usually around 5%, swell when they come into contact with water. Again, this is a purely physical aspect in which the metabolism of the seeds is not directly involved. For example, seeds which are already dead and no longer germinable swell just as well as those which are germinable.

Tasks



Swollen pea seeds



How much water do pea seeds absorb during swelling?

Investigate the increase in size of pea seeds due to swelling.



Equipment

Position	Material	Item No.	Quantity
1	Beaker, Borosilicate, tall form, 100 ml	46026-00	1
2	Graduated cylinder 100 ml, PP transparent	36629-01	1
3	Petri dish, d 100 mm	64705-00	1

Set-up

PHYWE excellence in science

excellence in science

- Fill a measuring cylinder of 100 ml capacity up to the 50 ml mark with tap water.
- Add 50 dry pea seeds and swirl the graduated cylinder back and forth several times to check for any air bubbles.
- Remove trapped air bubbles.
- $\circ~$ Read the level of the water in the graduated cylinder.



Remove any air bubbles

Procedure

- Pour the pea seeds with the water into a beaker of 100 ml capacity and cover it with half a Petri dish of 100 mm diameter.
- After two to three days, pour off the water from the beaker, take out the pea seeds and dry their surface carefully between an absorbent paper towel. The pea seeds have swelled.
- Fill the measuring cylinder again with tap water up to the 50 ml mark. Add the 50 swollen pea seeds and swirl the measuring cylinder back and forth several times to remove any trapped air bubbles. Read the water level.



Read water level





 $\mathbf{PH'}\mathbf{W}$

Report





Task 2	PHYWE xcellence in science
Choose the correct statements. Swelling is a purely physical process. Swelling is a purely chemical process. The water content of air-dry seeds is usually above 50%. That is why it is important that the see swell to increase their water content. The water content of air-dry seeds is usually around 5%. Check	eds
Task 3	PHYWE xcellence in science

Which statements are correct?			
The seed coat is permeable to water so that it can penetrate unhindered.			
Swelling can only take place under the influence of sunlight, as it derives its energy from photosynthesis.			
Through the hydration of the plasma in the cells of the embryo, a tugor pressure builds up, which is later the driving force for the germination of the seed.			
The seed coat is impermeable to water. It must be removed before swelling.			
Check			

PHYWE excellence in science

Task 4

Use your measurements to calculate the volume of the pea seeds before and after swelling.



Slide		Score/Total
Slide 13: Swelling		0/5
Slide 14: Swelling process		0/2
Slide 15: The seed coat		0/2
	Total	0/9
	Solutions	



www.phywe.de