


PHYWE Systeme GmbH & Co. KG  
Robert-Bosch-Breite 10  
D-37079 Göttingen

Telefon +49 (0) 551 604-0  
Fax +49 (0) 551 604-107  
E-mail info@phywe.com

### Operating instructions

 The unit complies with the corresponding EC guidelines







Fig. 1: 39843-01 Crystal lattice kit: sodium chloride (NaCl)

### Function and Applications

This kit includes all needed elements to build a Buckminster fullerene (C60) model.

Crystal lattice kit: sodium chloride NaCl					
Atoms	Holes (Angle)	Color	Dia. (mm)	Qty.	Drawing
Cl-Chlorine	6(90°)	Green	30	14	
Na-Sodium	6(90°)	Grey	23	13	
Links		Color	Length (mm)	Qty.	Drawing
Medium Bond		Grey white	44	54	
ABS plastic box (size: 23.5 × 17 × 3.5 cm)				1	